

## Handpiece | 手柄

NO CRACKING! 不开裂!



### Model:H1

连体式手柄  
频率自动跟踪 工作稳定  
Sealed handpiece  
Automatic frequency tracking,  
works steadily.



### Model:H2

接插式手柄 可高温高压消毒  
适用品牌: BAOLAI, EMS  
Detachable handpiece can be autoclaved  
under the high temperature and pressure.  
(135°C, 0.22MPa)  
Compatible with BAOLAI and EMS scalars



### Model:H3

铝合金材质 不开裂 更舒适  
Aluminum alloy material,  
No cracking! More comfortable!

接插式手柄 可高温高压消毒  
Detachable handpiece can be autoclaved  
under the high temperature and pressure.  
(135°C, 0.22MPa)

适用品牌: BAOLAI, EMS  
Compatible with BAOLAI and EMS scalars



## Handpiece | 手柄

宝莱铝合金手柄为全球首创,荣获15项国家专利,仿冒必究。  
Aluminum alloy handpiece is first invented by BaoLai  
in the world and awarded 15 national patents.  
Copy forbidden.




### 铝合金手柄

铝合金材质,反复高温高压消毒,  
手柄依然保持完整,不会开裂。  
有利于手柄散热,更好保护换能器,  
延长使用寿命。  
工作频率更稳定,为患者带来舒适感受。  
防滑设计,手握舒适,更高档次。

### Aluminum alloy handpiece

Made of aluminum, it can undergo repeating autoclave  
and remain uncracked.  
Aluminum better disperses heat, protects the transducer  
and therefore extends use life.  
Much more stable working frequency ensures patients' comfort.  
Anti-skidding design, more comfortable handling, higher grade



### 传统塑料手柄

经过反复高温高压消毒,手柄容易变黄、开裂。  
塑料材质不利散热,为患者带来酸痛感。  
散热不好带来换能器频率不稳定,影响使用寿命。

### Traditional plastic handpiece

Handpiece will turn brownish and be easy to crack after repeating autoclave.  
Plastic disperses heat very slow and therefore brings patients ache.  
The slow dispersing of heat also affects the frequency of transducer and use life.



## P5&P4&B5 Ultrasonic scaler | 超声洁牙机



### Model:P4

按钮式手柄 可高温高压消毒

Detachable handpiece can be autoclaved under the high temperature and pressure. (135°C, 0.22MPa)

主机输入 Rate: 24VAC 50/60Hz  
频率 Frequency: 28kHz±3kHz



### Model:B5

连体式手柄

频率自动跟踪 工作稳定

Sealed handpiece  
Automatic frequency tracking,  
works steadily.

主机输入 Rate: 24VAC 50/60Hz  
频率 Frequency: 28kHz±3kHz



### Model:P5

铝合金手柄 不开裂 更舒适

Aluminum alloy handpiece!  
No cracking! More comfortable!

按钮式手柄 可高温高压消毒

Detachable handpiece can be autoclaved under the high temperature and pressure. (135°C, 0.22MPa)

主机输入 Rate: 24VAC 50/60Hz  
频率 Frequency: 28kHz±3kHz



## P7&P6&B6 Ultrasonic scaler | 超声洁牙机

### SCALING/PERIO/ENDO

洁治功能 牙周功能 根管功能



### Model:P6

洁治功能 牙周功能 根管功能  
Scaling/Perio/Endo

按钮式手柄 可高温高压消毒

Detachable handpiece can be autoclaved under the high temperature and pressure. (135°C, 0.22MPa)

主机输入 Rate: 30VDC 50/60Hz  
频率 Frequency: 28kHz±3kHz



### Model:P7

洁治功能 牙周功能 根管功能  
Scaling/Perio/Endo

铝合金手柄 不开裂 更舒适

Aluminum alloy handpiece!  
No cracking! More comfortable!

按钮式手柄 可高温高压消毒

Detachable handpiece can be autoclaved under the high temperature and pressure. (135°C, 0.22MPa)

主机输入 Rate: 30VDC 50/60Hz  
频率 Frequency: 28kHz±3kHz



### Model:B6

洁治功能 牙周功能 根管功能  
Scaling/Perio/Endo

连体式手柄

频率自动跟踪 工作稳定

Sealed handpiece  
Automatic frequency tracking,  
works steadily.

主机输入 Rate: 30VDC 50/60Hz  
频率 Frequency: 28kHz±3kHz



# P9

## Auto-water Supply Scaler | 自动供水洁牙机 SCALING/ ENDO 洁治功能 根管功能

自动供水装置；  
电源内置；  
按钮式手柄，可高温高压消毒；  
铝合金手柄，荣获国家多项专利；  
洁治、根管功能；自动识别  
微电脑控制，频率自动跟踪。  
Auto-water supply device;  
Built-in power supply;  
Detachable handpiece can be autoclaved under  
the high temperature and pressure.(135°C,0.22MPa)  
Aluminum alloy handpiece;  
obtaining national patents;  
Scaling, Endo;Identifying scaler tip automatically;  
PC control, automatic frequency tracking.



Compact,space saving:243x145x150(mm)  
体积小，节约空间。

Water supply system without consumables, cost saving.  
供水系统零耗材，降低使用成本。

Aluminum alloy LED handpiece is available  
铝合金LED手柄可选配



H3 Aluminum alloy detachable handpiece  
H3铝合金按钮式手柄

Noiseless working  
静音工作

**NO Consumables!**  
**零耗材!**



Aluminum alloy LED handpiece is available  
铝合金LED手柄可选配



Functional shift buffer zone  
功能缓冲区



Prepositive water  
regulator  
水量开关装置



Filter design for  
protection  
过滤保护设计



Convenient for  
liquid adding  
方便添加液体



Built-in power supply  
电源内置



Convenient for  
liquid monitoring  
方便观察液体

# C7&C6&C5

Built-in scaler | 内置式洁牙机

Special design for dental unit 特别为牙椅设计



### Model:C6

按钮式手柄 可高温高压消毒

Detachable handpiece can be autoclaved under the high temperature and pressure. (135°C, 0.22MPa)

主机输入 Rate: 24VAC 50/60Hz  
频率 Frequency: 28kHz±3kHz



### Model:C5

连体式手柄 频率自动跟踪, 工作稳定。

Sealed handpiece Automatic frequency tracking, works steadily.

主机输入 Rate: 24VAC 50/60Hz  
频率 Frequency: 28kHz±3kHz



### Model:C7

铝合金手柄 不开裂 更舒适  
牙椅的最佳拍档

Aluminum alloy handpiece  
No cracking! More comfortable!  
The best match for dental unit

按钮式手柄 可高温高压消毒

Detachable handpiece can be autoclaved under the high temperature and pressure. (135°C, 0.22MPa)

主机输入 Rate: 24VAC 50/60Hz  
频率 Frequency: 28kHz±3kHz



# Wrench | 扳手

## Torque Wrench 限力扳手

Model:TW3

Aluminum alloy material.

It could hold the tips for at least 3000 times!

铝合金材质, 最少能夹持工作尖3000次!



Torque Wrench 限力扳手

Model:TW2

Plastic material  
塑料材质



Wrench 扳手

Model:TW1



Endo Wrench 根管扳手

Model:EW1

## Scaler tip (T series) | 工作尖 (T系列)



**T1**  
Length: 29.9MM  
N.W: 0.9g

去除龈缘上及齿间的结石。  
Removal of supragingival and interdental calculus.



**T2**  
Length: 28.6MM  
N.W: 0.9g

去除龈缘上的结石，尤其是难去除的牙结石及牙菌斑。  
Removal of heavy calculus and stains in the supragingival area.



**T3**  
Length: 28.6MM  
N.W: 0.9g

去除龈缘上及龈缘下的结石，尤其是去除狭小齿间的结石及牙周袋内的结石。  
Removal of supra and subgingival calculus. It provides easy access to interdental spaces and narrow pocket.



**T4**  
Length: 29.6MM  
N.W: 0.9g

去除龈缘上及齿间的结石。  
Removal of supragingival and interdental calculus.



**T5**  
Length: 31.8MM  
N.W: 1.0g

去除臼齿间的结石。  
Removal of calculus between molar teeth.



**T6**  
Length: 30.2MM  
N.W: 0.9g

去除龈缘、齿间及龈缘下的结石，此工作尖较长，尤其适用臼齿。  
Removal of supragingival and interdental calculus. This tip can be used in all quadrants and is very useful for the removal of hard calculus. Especially between molar teeth.



**T7**  
Length: 30.5MM  
N.W: 0.9g

去除龈缘上的较难去除的结石，此工作尖较长，尤其适用臼齿。  
Removal of supragingival calculus. This tip can be used in all quadrants and is very useful for the removal of hard calculus. Especially between molar teeth.



**T8**  
Length: 31.8MM  
N.W: 0.9g

适用于去除龈缘下及牙周袋内的结石。  
Removal of supra and subgingival calculus. It provides easy access to interdental spaces and narrow pocket.



**T9**  
Length: 32.5MM  
N.W: 1.0g

适用于去除龈缘下及牙周袋内的结石。  
Removal of supra and subgingival calculus. It provides easy access to interdental spaces and narrow pocket.



**E1**  
Length: 21.2MM  
N.W: 1.0g

⚠ E1 is compatible with BAOLAI, same as T series tips.

适用于有效地清洗和冲洗根管，不具备切削作用。  
Non-cutting E1 instrument for effective cleaning and irrigation of the root canal system.



宝莱洁牙机请使用宝莱品牌的工作尖。  
Please use Baolai's scaler tips for Baolai's scalers.

适用品牌 COMPATIBLE WITH:  
BAOLAI, EMS, MECTRON  
洁治Scaling: T1, T2, T4, T5, T6, T7  
牙周治疗Periodontics: T3, T8, T9  
根管Endo: E1



适用品牌  
COMPATIBLE WITH: SATELEC, NSK  
洁治Scaling: S1, S2, S4, S5, S6, S7  
牙周治疗Periodontics: S3, S8, S9



## Endo file | 根管锉

不锈钢材质，更适合根管荡洗。  
Stainless steel material,  
More suitable for cleaning and irrigation  
of the root canal system.  
长度 Length : 33MM

